

Traffi

FOR HANDS. FOR LIFE.

TD01 SUSTAIN

CARBON NEUTRAL BIODEGRADABLE NITRILE DISPOSABLE GLOVE

The TD01 is part of Traffi's Sustain disposable glove range. It's made from 100% premium quality nitrile with a 3.5g high stretch construction. It's powder free, food approved and has a textured finish for improved grip. It's also biodegradable, meaning you can dispose of uncontaminated gloves in the general waste stream with minimal planetary impact Key features

Key features

- Manufactured from 100% Nitrile
- Premium Quality 3.5g High stretch material
- Powder Free
- Textured finish for improved grip
- Made in Sri Lanka
- Food approved
- Carbon neutral and biodegradable

Sectors

- Healthcare
- Food Manufacturing & preparation
- Construction & utility
- Janitorial & cleaning
- Automotive
- Facilities maintenance

Standards:

- Tested to European Standards EN455: Medical Gloves for Single Use Parts 1,2,3&4
- Product Food Approved to EN1186
- EU Type Examination carried out by SATRA
- Tested to European Standards EN ISO 21420:2020 & EN ISO 374-1:2016+A1:2018 & EN ISO 374-5:2016



EN ISO 374-1: 2016/TYPE B



EN ISO 374-5: 2016



Dimensions

Size	Weight g/pc	Lgth (mm) Median	Palm width (mm)	Thickness		
				Cuff	Palm	Finger
Small	3.4 ± 0.2	→ 240	90±5	0.06 ± 0.01	0.08 ± 0.01	0.10 ± 0.01
Medium	3.5 ± 0.2	→ 240	100±5	0.06 ± 0.01	0.08 ± 0.01	0.10 ± 0.01
Large	4.2 ± 0.2	→ 240	110±5	0.06 ± 0.01	0.08 ± 0.01	0.10 ± 0.01
X Large	4.4 ± 0.2	→ 240	→ 120±5	0.06 ± 0.01	0.08 ± 0.01	0.10 ± 0.01



T +44 (0)1329311451



BIODEGRADES
BY **76.5%**
IN **539** DAYS



Medical



Pharmaceutical



Laboratory



General Purpose



Mechanical

Specifications

100 pieces per box

Size: XS- XXL

10 box per case

Colour: Blue

AQL: 1.5

Packaging

Size	Pack barcode	Case barcode
S	5055681046330	5055681046385
M	5055681046347	5055681046392
L	5055681046354	5055681046408
XL	5055681046361	5055681046415